


BDA063B

BLOOD COLLECTION SYSTEMS



TECHNICAL SPECIFICATIONS

| Bag(s) | Collection volume | Packaging | |
|------------|-------------------|-----------------------------|------------------|
| 2 PVC-DEHP | 450 mL +/- 10 % | 2 devices / peelable sachet | 80 devices / box |

Included items

Needle (16G)
Bactivam, Secuvam and Vacuvam
Break-Away Cannula in the Collection Line

Solution

63 mL CPDA1

Sterilisation

Steam

Labels

Universal pictograms (according to 3826-2)
ISBT128

Shelf life

24 months

Instructions for use

NOTIR01

Final products

Plasma
Red Cell Concentrate

Pishro Tashkhis Fardavar - Tehran - Iran

Tel : +9821-22262616-18 - Fax : +9821-26400055

www.fardavar.com - info@fardavar.org

Payam Special Economic Zone - Alborz - Iran

Reference: QA-US-GUI-06

Version: 00