


BQG263B

IN-LINE RED CELL FILTRATION



BLOOD PROCESSING



TECHNICAL SPECIFICATIONS

Bag(s)	Filter(s)	Collection volume	Packaging	
3 PVC-DEHP 1 PVC-TOTM	LEUCOFLEX LCR5	450 mL +/- 10 %	1 device / peelable sachet	12 devices / box

Included items

Needle (16G)
Bactivam, Secuvam and Vacuvam
Break-Away Cannula in the Collection Line

Solutions

63 mL CPD
100 mL SAGM

Sterilisation

Steam

Labels

Universal pictograms (according to ISO 3826-2)
ISBT128

Shelf life

2 years

Instructions for use

NOTIR02

Pishro Tashkhis Fardavar - Tehran - Iran

Tel : +9821-22262616-18 - Fax : +9821-26400055

www.fardavar.com - info@fardavar.org

Payam Special Economic Zone · Alborz · Iran

Final products

Plasma
Single Unit of Platelet Concentrate
Leucoreduced Red Cell Concentrate