


# BTA063B

## BLOOD COLLECTION SYSTEMS



### TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging	
2 PVC-DEHP 1 PVC-TOTM	450 mL +/- 10 %	2 device / peelable sachet	80 devices / box

#### Included items

Needle (16G)  
Bactivam, Secuvam and Vacuvam  
Break-Away Cannula in the Collection Line

#### Solutions

63 mL CPDA1

#### Sterilisation

Steam

#### Labels

Universal pictograms (according to ISO 3826-2)  
ISBT128

#### Shelf life

24 months

#### Instructions for use

NOTIR01

Pishro Tashkhis Fardavar - Tehran - Iran

Tel : +9821-22262616-18 - Fax : +9821-26400055

www.fardavar.com - info@fardavar.org


Payam Special Economic Zone - Alborz - Iran

#### Final products

Plasma  
Single Unit Platelet Concentrate  
Red Cell Concentrate

# BDA063B

## BLOOD COLLECTION SYSTEMS



### TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging	
2 PVC-DEHP	450 mL +/- 10 %	2 devices / peelable sachet	80 devices / box

#### Included items

Needle (16G)  
Bactivam, Secuvam and Vacuvam  
Break-Away Cannula in the Collection Line

#### Solution

63 mL CPDA1

#### Sterilisation

Steam

#### Labels

Universal pictograms (according to 3826-2)  
ISBT128

#### Shelf life

24 months

#### Instructions for use

NOTIR01

#### Final products

Plasma  
Red Cell Concentrate

Pishro Tashkhis Fardavar - Tehran - Iran

Tel : +9821-22262616-18 - Fax : +9821-26400055

www.fardavar.com - info@fardavar.org


Payam Special Economic Zone - Alborz - Iran

Reference: QA-US-GUI-06

Version: 00

# BSA063B

## BLOOD COLLECTION SYSTEMS



### TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging	
1 PVC-DEHP	450 mL +/- 10 %	1 device / peelable sachet	90 devices / box

#### Included items

Needle (16G)  
Bactivam, Secuvam and Vacuvam  
Break-Away Cannula in the Collection Line

#### Solution

63 mL CPDA1

#### Sterilisation

Steam

#### Labels

Universal pictograms (according to ISO 3826-2)  
ISBT128

#### Shelf life

18 months

#### Instructions for use

NOTIR01

#### Final product

Whole Blood

Pishro Tashkhis Fardavar - Tehran - Iran

Tel : +9821-22262616-18 - Fax : +9821-26400055

www.fardavar.com - info@fardavar.org

Payam Special Economic Zone - Alborz - Iran

Reference: QA-US-GUI-07


Version: 00

# BQG263B

## IN-LINE RED CELL FILTRATION



BLOOD PROCESSING



### TECHNICAL SPECIFICATIONS

Bag(s)	Filter(s)	Collection volume	Packaging	
3 PVC-DEHP 1 PVC-TOTM	LEUCOFLEX LCR5	450 mL +/- 10 %	1 device / peelable sachet	12 devices / box

#### Included items

Needle (16G)  
Bactivam, Secuvam and Vacuvam  
Break-Away Cannula in the Collection Line

#### Solutions

63 mL CPD  
100 mL SAGM

#### Sterilisation

Steam

#### Labels

Universal pictograms (according to ISO 3826-2)  
ISBT128

#### Shelf life

24 months

#### Instructions for use

NOTIR02

Pishro Tashkhis Fardavar - Tehran - Iran

Tel : +9821-22262616-18 - Fax : +9821-26400055

www.fardavar.com - info@fardavar.org

Payam Special Economic Zone - Alborz - Iran

#### Final products

Plasma  
Single Unit of Platelet Concentrate  
Leucoreduced Red Cell Concentrate