


BSD021B

ACDA Solution



TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging
1 PVC-DEHP	450 mL +/- 10%	1 device / peelable sachet 18 (Small) – 72 (Big) devices/box

Included Items

- Junk (CJ0325G)

Final Products

- Whole Blood

Solutions

- 250 mL ACDA

Sterilization

- Steam

Labels

- Universal pictograms (according to ISO 3826-1)
- ISBT128

Shelf Life

- 24 Months


Instructions for Use

- NOTIRO7

Pishro Tashkhis Fardavar - Tehran-Iran
 Tel: +9821-22262616-18 - Fax: +9821-26400055
 www.fardavar.com - info@fardavar.org
 Payam Special Economic Zone Alborz Iran

BSC021B

Sodium Citrate 4% Solution



TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging
1 PVC-DEHP	450 mL +/- 10%	1 device / peelable sachet 18 (Small) – 72 (Big) devices/box

Included Items

- Junk (CJ0325G)

Final Products

- Whole Blood

Solutions

- 250 mL Sodium Citrate 4%

Sterilization

- Steam

Labels

- Universal pictograms (according to ISO 3826-1)
- ISBT128

Shelf Life

- 24 Months

Instructions for Use


- NOTIRO5

Pishro Tashkhis Fardavar - Tehran-Iran
 Tel: +9821-22262616-18 - Fax: +9821-26400055
 www.fardavar.com - info@fardavar.org
 🏢 Payam Special Economic Zone Alborz Iran

LSA062A

BLOOD COLLECTION SYSTEMS

Letting Blood Bag



TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging
1 PVC-DEHP	450 mL +/- 10%	1 device / peelable sachet 25 (Small) – 90 (Big) devices/box

Included Items

- Needle (G16)
- Break-Away Cannula in the Collection Line

Solutions

- 35 mL CPDA1

Sterilization

- Steam

Labels

- Universal pictograms (according to ISO 3826-1)
- ISBT128

Shelf Life

- 24 Months

Instructions for Use

- NOTIRO3

Final Products

- Whole Blood