


# LSA062A

## BLOOD COLLECTION SYSTEMS

### Letting Blood Bag



## TECHNICAL SPECIFICATIONS

Bag(s)	Collection volume	Packaging
1 PVC-DEHP	450 mL +/- 10%	1 device / peelable sachet 25 (Small) – 90 (Big) devices/box

#### Included Items

- Needle (G16)
- Break-Away Cannula in the Collection Line

#### Solutions

- 35 mL CPDA1

#### Sterilization

- Steam

#### Labels

- Universal pictograms (according to ISO 3826-1)
- ISBT128

#### Shelf Life

- 24 Months

#### Instructions for Use

- NOTIRO3

#### Final Products

- Whole Blood